

TranspoNet IntelliSend

The IntelliSend Station protects high-value items like blood bags and medication with a secure safety drawer, ensuring that only **authorized personnel** can access them—ideal for wards with frequent visitor traffic.



Station Setup

1 IntelliSend Station

Optional without:

2 Drive tube connection


3 Beveled cover plate

4 Lockable base cabinet with drawer

5 Upper cover with service flap

6 SafetyLock

7 Standard collection basket



Features:

- » Easy integration into niches or tight spaces (60 cm wide)
- » Secured storage for up to six carriers
- » Illuminated door indicating system & station status
- » Energy-efficient components and automatic power-saving mode
- » Card reader with multi-factor authentication
- » RFID-enabled for tracking and empty carrier management
- » Modular system allows on-site upgrades
- » Ergonomic display
- » Use as a pass-through or end station

Customer Benefits

- » Ergonomic handling reduces work-related strain
- » Clear visual status indicators improve usability
- » Secure and traceable sending/receiving via RFID
- » Multi-factor authentication for enhanced security
- » Easy maintenance with front-accessible components
- » Optional mobile notifications for incoming items
- » Flexible installation options to suit various environments
- » Safe operation with pinch protection and motorized door

CONTACT



Swisslog Healthcare
Webereiweg 3
5033 Buchs
Switzerland
T: +41 62 837 41 41
contact.ch@swisslog-healthcare.com

